

LAP-5® Disposable Endoscopic Instruments
Instructions for use

Ref. no.: 0208-DSLAP501, 0208-DSLAP502, 0208-DSLAP503, 0208-DDLAP501, 0208-DGLAP501, 0208-DGLAP502, 0208-DGLAP503, 0208-DGLAP504, 0208-DGLAP505, 0208-DGLAP501B, 0208-DSLAP502B, 0208-DSLAP503B, 0208-DDLAP501B, 0208-DGLAP501B, 0208-DGLAP502B, 0208-DGLAP503B, 0208-DGLAP504B, 0208-DGLAP505B

 <p>KONMEX Konmex Sp. z o.o. ul. Bartycka 63A/6, 00-716 Warszawa, Poland</p>	<p>Contact information: Phone: +48 22 730 13 94 Fax: +48 22 780 06 09</p>		<p align="center">ENG KX-IFU-017-ENG_03</p>
---	---	--	--

Attention:

Read and understand all information contained in these instructions for use. Failure to do so properly may lead to serious surgical consequences. This instructions for use cannot be used as a manual for surgical techniques used in minimal invasive surgery. To learn adequate knowledge about surgical technique it is necessary to contact our company or authorized distributor and to acquaint with appropriate technical instructions, professional medical literature and graduate proper training under supervision of surgeon experienced in techniques of endoscopic surgery.

Instrument is delivered sterile and is intended for single use.

Indications:

Disposable endoscopic instruments are indicated for cutting, grasping, dissecting and coagulation of tissue in laparoscopic and thoracoscopic surgical procedures. They are intended for single patient and procedure use.

Intended users: Disposable Endoscopic Instruments are restricted to use exclusively by qualified medical staff.

Patient target group: Target age groups include all the adult and young patients, males and females.

Contraindications:


The use of disposable endoscopic instruments is contraindicated whenever endoscopic surgical techniques are contraindicated for any reason.

Prior to the use:

Carefully inspect the shipping carton, its contents and individual pouch for any sign of damage. If damage is visible, do not use the instrument.

Illustration of the instrument (pic. I):

- A. Jaws
- B. Shaft
- C. Rotating knob
- D. HF connector
- E. Thumb handle
- F. Front handle
- G. Ratchet trigger
- H. Ratchet bar




Pic. I

Pic. II

Instructions for use:

1. Open the package by using standard aseptic technique.
2. Ensure that the product functions and is intact.
3. Remove paper protectors from jaws and rotating knob.
4. If ratcheted instrument is used, open jaws and push trigger backwards in order to engage the ratchet mechanism (pic. II).
5. Close the handles to the desired grasping position. The instrument remains locked on the tissue (pic. II).
6. Move the trigger forward to release the jaws. (pic. III).
7. In order to use ratcheted instrument as a non-ratcheted one leave the trigger in the forward position. The instrument will open and close freely (pic. III).
8. Use rotating knob to turn instrument jaws in any direction (pic. IV)



Pic. III

Pic. IV

Electrosurgery:

First, connect the electrosurgical cord (not furnished with the instrument) to the device by placing the female end of the cord on the 4 mm male HF connector. Plug the other end of the cord into the monopolar receptacle of the HF generator. Attach return electrode to the patient body and connect it with the respective receptacle of the HF generator. If instrument and/or return electrode is not properly connected to the generator electrosurgery will not be possible to perform. Recommended maximum output power of the generator to be used with the device is 350W for cut and 120W for coagulation with blend cut power between above values.

Rated accessory voltage of the device – 1 500V.

⚠️ Electrocauter precautions:

1. A complete understanding of the principle of monopolar electrocauter procedures is necessary to avoid accidental shocks, burns, or potential gas embolism to the patient.
2. Be sure that the entire area of the return electrode has been properly attached to the patient's body and is as close to the operating field as possible. Incomplete body-electrode contact may lead to burns and/or inability to perform electrocauter.
3. The patient should not come into contact with metal parts which are earthed or which have an appreciable capacitance to earth (for example operating table supports, etc.) as it may lead to burn injury of the patient. The use of antistatic sheeting is recommended for patient protection.
4. To protect patient from burns skin-to-skin contact (for example between the arms and body of the patient) should be avoided, for example by insertion of dry gauze.
5. The use of flammable anaesthetics or oxidizing gases such as nitrous oxide (N₂O) and oxygen should be avoided if a surgical procedure is carried out in the region of the thorax or the head, unless these agents are sucked away. Combustible gases can ignite during electrocauter, seriously injuring the patient and the surgeon.
6. Non-flammable agents should be used for cleaning and disinfection wherever possible. Flammable agents used for cleaning or disinfecting, or as solvents of adhesives, should be allowed to evaporate before the application of HF surgery. There is a risk of pooling of flammable solutions under the patient or in body depressions such as the umbilicus, and in body cavities such as the vagina. Any fluid pooled in these areas should be mopped up before HF surgical instrument is used. Residual flammable agents can ignite during HF surgery, leading to severe thermal injuries of the patient and the surgeon.
7. Attention should be called to the danger of ignition of endogenous gases. Some materials, for example cotton, wool and gauze, when saturated with oxygen may be ignited by sparks produced in normal of the HF surgical instrument leading to thermal injuries of the patient and the surgeon.
8. For patients with cardiac pacemakers or other active implants, a possible hazard exists because interference with the action of the pacemaker may occur, or the pacemaker may be damaged. In case of doubt, approved qualified advice should be obtained.
9. If any physiological monitoring equipment is used simultaneously with HF generator on the same patient, any monitoring electrodes (including monitoring device) should be placed as far as possible from HF generator. Needle monitoring electrodes are not recommended as they may cause patient burns. The use of monitoring systems incorporating high frequency current limiting devices is recommended.
10. The cables to the electrocauter instruments (including HF generator) should be positioned in such a way that contact with the patient or other leads is avoided to prevent short-circuit or patient burns in case of insulation damage.
11. Temporarily unused electrocauter instruments (including HF generator) should be stored in a location that is isolated from the patient.
12. For surgical procedures where the HF current could flow through parts of the body having a relatively small cross sectional area, the use of bipolar or pure heat techniques may be desirable in order to avoid unwanted coagulation.
13. Do not activate the generator until instruments' jaws are in contact with tissue or are in a position to deliver high frequency energy to the tissue. Premature activation can lead to coagulation at unintended sites.
14. Keep the output power as low as possible to achieve the desired effect. Surgeon is fully responsible for the correct coagulation time and power. Prolonged coagulation time and/or excessive power may lead to tissue charring and widening of the area of lateral lesions.
15. Avoid HF output settings of the generator where maximum output voltage may exceed rated accessory voltage. Exceeding the rated voltage may damage the insulation and result in thermal injury of the patient and the operator.
16. Apparent low output or failure of the HF surgical equipment to function correctly at the normal operating settings may indicate faulty application of the neutral electrode or poor contact in its connections. In this case, the application of the neutral electrode and its connections should be checked before selecting a higher output power.
17. When using electrocauter, verify that jaws of the instrument are not in contact with a conductive irrigation fluid. HF current flowing through conductive fluid may lead to burns in multiple areas inside the patient body.
18. Electrocauter generators used with these devices may cause unintended destruction of tissue and are dangerous if operated improperly. Read carefully instructions for use of the generator prior to the procedure.
19. Sufficient care and distance must be maintained during use to prevent arcing to other instruments leading to unintended coagulation of the sites remaining in direct contact with these instruments.

⚠️ Additional warnings and precautions:

1. Minimally invasive surgery should be performed only by physicians who are thoroughly trained in minimally invasive techniques.
2. To avoid injury to internal organs, a pneumoperitoneum must be maintained during the use of disposable endoscopic instrumentation.
3. Verify that the devices are compatible with other products that will be used in surgery prior to the procedure. Incompatibility may lead to an extended procedure time, inability to perform surgery or necessity to convert to an open surgery.
4. Dispose of all opened instruments whether used or unused to prevent accidental use of a contaminated device.
5. Use immediately after opening. Keeping the instruments post package opening leads to their contamination and creates a risk of an infection to the patient.
6. This product is intended for single patient and procedure use. Resterilization, reuse, modification may lead to serious consequences with death of patient included.
7. Take care to discard the product and packing after use, as well as unused but opened devices in accordance with hospital waste disposal practices and local regulations including, without limitation, those pertaining to human health and safety and the environment.
8. If any serious incident that has occurred in relation to the device it should be reported to the manufacturer and the competent authority of the Member State, in which the user and/or patient is established.



Keep dry



Consult electronic instructions for use



Manufacturer



Date of manufacture



Caution



Do not resterilize



Do not use if package is damaged and consult instructions for use



Use-by date



Sterilized using ethylene oxide



Catalogue number



Batch code



Quantity in package



Single sterile barrier system



Do not re-use



Medical Device



Unique device identifier

*The hard copies of instructions for use delivered with Konmex products are always in english language.
If you require a hard copy of IFU in other language, you can contact Konmex Sp. z o. o.
at regulatory@konmex.com or **+48 (22) 730 13 94**.*

*Please scan the below QR code with the appropriate application.
It will connect you with Grena Ltd. website where you can choose eIFU in your preferable language.*

You can enter the website directly by typing in www.grena.co.uk/IFU in your browser.

*Make sure that paper version of IFU in your possession is in the latest revision prior to use of the device.
Always use the IFU in the latest revision.*

