## NE'X Glue<sup>®</sup> Surgical Adhesive

### Information for patients regarding NE'X Glue<sup>®</sup> Surgical Adhesive

#### General description of the device

NE'X Glue<sup>®</sup> Surgical Adhesive is a two-component surgical glue composed of purified bovine serum albumin (BSA) solution and glutaraldehyde solution. Two main components are separated in double-chambered syringe. The solutions are dispensed and mixed by applicator tips, which are designed to provide precise and reproducible mixing of components during application.

NE'X Glue<sup>®</sup> Surgical Adhesive is indicated to bond, seal and/or reinforce soft tissue. Can be applied as an adjunct to staples, sutures, electrocautery or patches as well as alone for sealing or reinforcing parenchymal organs when other standard methods are impractical or ineffective. Another application is surgical mesh fixation in hernia surgery. Soft tissues where NE'X Glue<sup>®</sup> is effective are vascular, cardiac, pulmonary, dural, esophageal, gastric, intestinal, colorectal, pancreatic, splenic, biliary, hepatic and genitourinary.

#### Materials and substances

NE'X Glue<sup>®</sup> Surgical Adhesive is available in 3 volumes – 2mL, 5mL, and 10mL. Each syringe is composed of bovine serum albumin (BSA) and glutaraldehyde solutions in a 4:1 ratio.

The specification for the BSA solution is a 45% (weight/volume ratio) solution. The amount of raw BSA needed for each size are: 1,04g (2ml syringe), 1,5g (5ml syringe), and 3,71g (10ml syringe).

The specification for the glutaraldehyde solution is a 10% (weight/volume ratio) solution. The maximum quantity of glutaraldehyde coming into contact with the patient when using a single device is 0,58ml (2ml syringe), 1,03ml (5ml syringe), and 2,06ml (10ml syringe) for each configuration.

#### Warnings, precautions or actions to be taken by the Patient or healthcare professional in the possibility of being influenced by external factors, medical examinations or environmental conditions

The warnings and precautions for healthcare professionals are described in the Instructions for USE (IFU) available at [www.grena.co.uk/IFU](http://www.grena.co.uk/IFU). NE'X Glue<sup>®</sup> Surgical Adhesive should not be used in case of known sensitivity to materials of bovine origin. Possible side effects may include, but are not limited to: failure of glue to adhere to tissue, inflammatory response, immune response, allergic reaction, application to tissue not targeted for the procedure, tissue necrosis, vessel obstruction, bronchus obstruction, luminal obstruction, tissue mineralization, thrombosis and thromboembolism, pulmonary emboli, injury to vessels or tissue, transmission of infectious agents of animal origin.

NE'X Glue<sup>®</sup> Surgical Adhesive is **safe** for diagnostic procedures and imaging using Magnetic Resonance (**MR**) **technology** (i.e., an item that poses no hazard in all MR environments).

#### Expected lifetime

NE'X Glue<sup>®</sup> Surgical Adhesive is a long- term implantable device, stays in the body for more than 30 days. Device is degraded via proteolysis (breakdown of proteins in the body). The rate of degradation will depend on the location and the amount of product applied. NE'X Glue<sup>®</sup> Surgical Adhesive degrades more quickly in areas with more blood vessels and less quickly in areas where there are less blood vessels. Thin applications of NE'X Glue<sup>®</sup> Surgical Adhesive will degrade more quickly than thick applications.

**Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.**

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